


Moist Heat Therapy Water Temperature Testing and Cleaning Log


Please allow adequate time for water to reach working temperature prior to testing.

If thermostat must be adjusted or replaced, remove unit from operation. Temperature must be re-tested prior to use of hot packs.

	Date	Time	Temp.	Corrective Action **	Drained & Cleaned?	Staff Initials
Mon						
Tues						
Wed						
Thurs						
Fri						
Sat						
Sun						
Mon						
Tues						
Wed						
Thurs						
Fri						
Sat						
Sun						
					 Unit cleaned?	



Temperature
 should not
 exceed
71-74°C*

	Date	Time	Temp.	Corrective Action **	Drained & Cleaned?	Staff Initials
					 Unit cleaned?	

****Corrective Action**
 Is the temperature within appropriate range? If not:
 "A" = Thermostat adjusted and re-tested until specific temperature is reached.
 "R" = Thermostat Replaced

*** Manufacturer Recommended**


Moist Heat Therapy Water Temperature Testing and Cleaning Log

Please allow adequate time for water to reach working temperature prior to testing.

If thermostat must be adjusted or replaced, remove unit from operation. Temperature must be re-tested prior to use of hot packs.

	Date	Time	Temp.	Corrective Action **	Drained & Cleaned?	Staff Initials
Mon						
Tues						
Wed						
Thurs						
Fri						
Sat						
Sun						
Mon						
Tues						
Wed						
Thurs						
Fri						
Sat						
Sun						
					 Unit cleaned?	



	Date	Time	Temp.	Corrective Action **	Drained & Cleaned?	Staff Initials
					 Unit cleaned?	

****Corrective Action**

Is the temperature within appropriate range? If not:
 "A" = Thermostat adjusted and re-tested until specific temperature is reached.
 "R" = Thermostat Replaced

*** Manufacturer Recommended**